

ILS Laboratories

8222 Vickers St, Suite 106, San Diego, CA 92111
(619) 329-3999 | ils-lab.com

Selank 10mg - 10mg

Tested for: Felix Chem LLC

PASS

COA #: COA-2026-EP_0-C
Lot Number: 052826-SLK-1
Accession #: ACC-2026-3606
Labeled Content: 10mg

Method:
Analysis Date: 06/07/2026
Appearance: Good
Sample Matrix: Lyophilized
Date Received: 05/30/2026



Scan to verify
authenticity at ils-lab.com
Access Code: F2P2NKFZ

Identity

Selank 10mg

Peptide Purity

99.73%



Fentanyl
Free


Full QC Panel

Analyte	Specification	Result	Unit	Status
Peptide Purity (HPLC)	>= 95.0%	99.73%	%	PASS
Net Peptide Content	Report Only	10.25	mg	N/A
Identity (HPLC-RTM)	Selank	Confirmed	-	PASS
Heavy Metals				
Arsenic (As)	<= 1.5 ppm	0.224 ppm	ppm	PASS
Cadmium (Cd)	<= 0.5 ppm	Not Detected	ppm	PASS
Lead (Pb)	<= 1.0 ppm	0.054 ppm	ppm	PASS
Mercury (Hg)	<= 1.5 ppm	Not Detected	ppm	PASS
Chromium (Cr)	<= 10.0 ppm	0.21 ppm	ppm	PASS
Fentanyl Screen	Immunoassay, 50 ng/mL cutoff	Not Detected	-	PASS



Selank 10mg 10mg - 052826-SLK-1

HPLC Chromatogram



Selank 10mg 10mg - 052826-SLK-1: UV Chromatogram

Sterility Testing (PCR)

Test	Specification	Result	Status
Sterility (PCR)	No Growth	No Growth	PASS



Handwritten signature

Dr. Greg Kalyuzhny
Lab Director
6/7/2026

COA #: COA-2026-EP_0-C
Access Code: F2P2NKFZ
Verify: portal.ils-lab.com/verify/S0quA6NWPzltkDf
Issued: 6/7/2026

Endotoxin Testing (USP <85>)

Test	Specification	Result	Status
Endotoxin (USP <85>)	Report Result	NMT 0.05 EU/mL	Reported
About this result: Endotoxin is reported as a quantitative value. Acceptable limits vary by product type and matrix, so no universal pass/fail threshold applies to RUO products. This result is below commonly referenced endotoxin thresholds.			

Notes & Methodology

- 1. Date Tested: 06/07/2026. Methods: Full QC Panel.
- 2. The sample was confirmed to be Selank 10mg by HPLC. Identification by chromatographic retention time comparison with a reference standard.
- 3. Endotoxin tested per USP <85> kinetic turbidimetric method. Acceptance criteria per client specification.
- 4. Peptide purity determined by RP-HPLC area normalization at 214 nm. Value represents the percentage of the target peptide relative to all peptide-related peaks. Non-peptide process-related impurities, if detected, are excluded from the calculation.




Dr. Greg Kalyuzhny
Lab Director
6/7/2026

COA #: **COA-2026-EP_0-C**
Access Code: **F2P2NKFZ**
Verify: <portal.ils-lab.com/verify/S0quA6NWPzltkDf>
Issued: 6/7/2026